



## Gloreha (Hand Rehabilitation Glove): summary of clinical results

August 2013

Although Gloreha (Hand Rehabilitation Glove) is a new device, its effectiveness has already been verified thanks to some clinical trials which involved:

- **post-stroke patients in the sub-acute phase** (up to a maximum of 30 days after the event;  $MAS \leq 3$ )<sup>1,2</sup>;
- **post-stroke patients in the chronic phase** (at least 3 months after the event;  $MAS \geq 1+$ ;  $MAS \leq 3$ ) treated with botulinum toxin<sup>3,4,2</sup>;
- **post-stroke patients with left neglect or general attention deficit** (cerebral stroke of the right hemisphere;  $MAS < 3$ )<sup>5,2</sup>;

---

<sup>1</sup> Vanoglio F., Luisa A., Garofali F., Mora C.; *Evaluation of the effectiveness of Gloreha (Hand Rehabilitation Glove) on hemiplegic patients. Pilot study*; presented at XIII Congress of Italian Society of Neurorehabilitation, 18-20 April 2013, Bari (Italy).

<sup>2</sup> Parrinello I., Faletti S., Santus G.; *Use of a continuous passive motion device for hand rehabilitation: clinical trial on neurological patients*; submitted to 41° National Congress of Italian Society of Medicine and Physical Rehabilitation, 14-16 October 2013, Rome (Italy).

<sup>3</sup> Stagno D., Baricich A., Invernizzi M., Grana E., Cisari C.; *Use of a robotic device in the rehabilitation treatment after botulinum toxin (type A) injection on spastic upper limb after stroke. Pilot study*; presented at XIII Congress of Italian Society of Neurorehabilitation, 18-20 April 2013, Bari (Italy).


<sup>4</sup> Molteni F., Mulè C., Caimmi M., Taveggia G., Missud S., Bianchi P.; *Use of continuous passive motion device for the hand on post-stroke patients: pilot study*. Presented at 39° National Congress of Italian Society of Medicine and Physical Rehabilitation, 16-19 October 2011, Reggio Calabria (Italy).

<sup>5</sup> Varalta V., Smania N., Geroin C., Fonte C., Gandolfi M., Picelli A., Munari D., Ianes P., Montemezzi G., La Marchina E.; *Effects of passive rehabilitation of the upper limb with robotic device Gloreha on visual-spatial and attentive exploration capacities of patients with stroke issues*; presented at Congress "Riabilitazione: una scienza in cammino", 18-20 March 2013, La Villa (Bolzano - Italy) and XIII Congress of Italian Society of Neurorehabilitation, 18-20 April 2013, Bari (Italy).


- **other neurological patients:** patients affected by multi-neurite, traumatic brain injury, meningioma of the spinal cord<sup>2</sup>.

Here are the main conclusions of the clinical trials on Gloreha:

- Gloreha can **improve functional independence** of post-stroke patients in the sub-acute phase.




- Gloreha can **increase grip and pinch strength** on the paretic side of neurological patients in the sub-acute phase.




- Gloreha can improve unimanual and bimanual coordination and dexterity of neurological patients, especially in the sub-acute phase.

**Vanoglio et al., 2013.** The evaluation of unimanual (paretic and healthy limb) and bimanual coordination and dexterity shows averages of delta t1-t0 (t1: after the treatment; t0: before the treatment) in the group treated with Gloreha higher than in the control group. The patients were post-stroke in the sub-acute phase. The treatment consisted of 30 sessions; the duration of each session was 40 minutes. The patients included in the control group were treated by therapists for the same time as the group treated with Gloreha.




**Varalta et al., 2013.** At the pre-treatment evaluation, after 4 months since hemorrhagic stroke, the patient F.M. had high difficulty in using the left hand to perform exercises that require good capacities of manual dexterity. After the treatment with Gloreha, the patient showed a small improvement in both tests (Purdue PegBoard Test and Nine Hole Peg Test). The treatment was made up of 10 rehabilitation sessions of 25 minutes.




**Parrinello et al., 2013.** After the treatment with Gloreha, 4 out of 12 neurological patients showed an improvement of manual ability and dexterity, evaluated by Nine Hole Peg Test and Abilhand. The treatment with Gloreha consisted of 10 daily sessions (30 minutes per session); the patients started the treatment at least 3 weeks since the event.

- Gloreha can **increase the motricity (active movement)** of neurological patients.



- Gloreha can **decrease or prevent the hypertonia** of neurological patients.



- Gloreha can **reduce the oedema** of neurological patients.


**Parrinello et al., 2013.** After the treatment with Gloreha two neurological patients showed a reduction of oedema. The treatment with Gloreha consisted of 10 daily sessions (30 minutes per session); the patients started the treatment at least 3 weeks since the event.

**Vanoglio et al., 2013.** All the patients, treated with Gloreha or by a therapist, decreased the level of oedema. The treatment consisted of 30 sessions; the duration of each session was 40 minutes.

**Molteni et al., 2011.** After the treatment with the prototype of Gloreha 5 out of 7 post-stroke patients with oedema showed a reduction of their symptoms.

- Gloreha can **improve the capacities of visual-spatial exploration** in post-stroke patients with neglect.


**Varalta et al., 2013.** After the treatment with Gloreha, the patient F.M. showed an improvement in the Line Crossing Test, Bells Test, Line Bisection Test.





**Parrinello et al., 2013.** After the treatment with Gloreha, the patients with neglect showed a partial improvement in their capacity of visual-spatial exploration. The treatment with Gloreha consisted of 10 daily sessions (30 minutes per session); the patients started the treatment at least 3 weeks since the event.

- Gloreha can **increase the speed of noticing visual stimuli in the left hemifield and reduce the reaction time during sustained attention tests** in post-stroke patients with neglect or general attention deficit.

**Varalta et al., 2013.** All the patients involved in the clinical test showed a decrease of the reaction time in the stimuli test (Saccadic Training).




Two patients out of four showed an improvement in SART (Sustained Attention to Response Task) Test, both in the reaction time and the score. The treatment with Gloreha was composed of 10 sessions (25 minutes per session).



- Gloreha is **well tolerated by patients**: all the tests show a high compliance level.

- **Gloreha can increase the time of treatment, at a limited cost:** Gloreha can bring a significant advantage in economic terms and enable the reallocation of human resources (i.e. therapists).



Vanoglio et al., 2013 concludes that **Gloreha treatment is not inferior to the standard treatment.** Actually the group treated with Gloreha showed in general more significant clinical improvements in comparison with the control group, but further trials are needed to have more solid data.

New clinical trials are in progress to substantiate and expand the results mentioned above, using a bigger sample size and also testing Gloreha on patients with:

- multiple sclerosis;
- brain tumour;
- spinal cord injury;
- traumatic brain injury.

Gloreha was developed in accordance with the conclusions of the most recent scientific literature:

- **a rehabilitation treatment after stroke has to be intensive, repetitive, functional, task-oriented, specific, customizable<sup>6,7</sup>;**

<sup>6</sup> Kollen BJ, Lindeman E. *Understanding the pattern of functional recovery after stroke: Facts and theories.* Restorative Neurology and Neuroscience 2004;22:281–299

<sup>7</sup> Van Peppen RP, Kwakkel G, Wood-Dauphinee S, Hendriks HJ, Van der Wees PJ, Dekker J. *The impact of physical therapy on functional outcomes after stroke: what's the evidence?* Clin Rehabil 2004;18:833–62.

- biomechatronic devices, robotic systems and computer interfaces can be useful for facilitating functional recovery and reducing motor impairment of upper limb<sup>8,9,10,11</sup>;
- a repetitive training program based on flexion and extension of the fingers can improve the function of plegic hand after stroke: it can reduce the oedema, maintain or increase the range of motion and muscle lengths, prevent spasticity<sup>12,13,14,15,16</sup>;
- passive stretch can prevent connective tissue accumulation in muscle<sup>17</sup> and help motor learning<sup>18,19,20</sup>;
- the exercises with devices works on plastic processes in the central nervous system thanks to neuromotor, visual and audio feedback: **multi-sensory action-observation systems enable patients to re-learn impaired motor functions through the activation of internal action-related representations**<sup>21,22</sup>;

---

<sup>8</sup> Masiero S, Celia A, Rosati G, Armeni M. *Robotic-assisted rehabilitation of the upper limb after acute stroke*. Arch Phys Med Rehabil 2007;88:142-9.

<sup>9</sup> Kiguchi K, Iwami K, Yasuda M, Watanabe K, Fukuda T. *An exoskeleton robot for human shoulder joint motion assist*. IEEE/ASME Trans Mech 2003;8:126-36.

<sup>10</sup> Hesse S, Schulte-Tiggens G, Konrad M, Bardeleben A, Werner C. *Robot-assisted arm trainer for the passive and active practice of bilateral forearm and wrist movement in hemiparetic subjects*. Arch Phys Med Rehabil 2003;84:915-20.

<sup>11</sup> Volpe B, Krebs H, Hogan N, Edelstein L, Diels C, Aisen M. *A novel approach to stroke rehabilitation: Robot aided sensorimotor stimulation*. Neurology 1999;54:1938-44.

<sup>12</sup> Carey JR, Kimberley TJ, Lewis SM, et al. *Analysis of fMRI and finger tracking training in subjects with chronic stroke*. Brain 2002;125(pt 4):773-88.

<sup>13</sup> Carey JR, Durfee WK, Bhatt E, et al. *Comparison of finger tracking versus simple movement training via telerehabilitation to alter hand function and cortical reorganization after stroke*. Neurorehabil Neural Repair 2007;21:216-32.

<sup>14</sup> Giudice M.L. *Effects of Continuous Passive Motion and elevation on hand edema*. The American Journal of Occupational Therapy, 1990; Oct; 44(10):914-21

<sup>15</sup> Dirette D., Hinojosa J. *Effects of continuous passive motion on the edematous hands of two persons with flaccid hemiplegia*. The American Journal of Occupational Therapy, 1994; May; 48(5):403-409.

<sup>16</sup> Xu X.L., Tong R.K., Song R., Zheng X.J., Leung W. W. *A comparison between electromyography-driven robot and passive motion device on wrist rehabilitation for chronic stroke*. Neurorehabilitation and neural repair, 2009; Jun 16

<sup>17</sup> Williams P E, Catanese T, Lucey E G, Goldspink G: *The importance of stretch and contractile activity in the prevention of connective tissue accumulation in muscle*. J Anat 1988; 158: 109-114.

<sup>18</sup> Seitz et al. *A common frame work for perceptual learning*. Curr. Neurobiol 2007, 17:1-6.

<sup>19</sup> Winter et al. *Hands-on therapy intervention for upper limb motor dysfunction following stroke*. Cochrane Database Syst Rev 2011 Jun 15; 6:CD006609.

<sup>20</sup> Hwang CH, Seong JW, Son DS. *Individual finger synchronized robot-assisted hand rehabilitation in subacute to chronic stroke: a prospective randomized clinical trial of efficacy*. Clinical Rehabilitation 2012;26(8); 696–704.

<sup>21</sup> Lee MM, Cho HY, Song CH. *The mirror therapy program enhances upper-limb motor recovery and motor function in acute stroke patients*. Am J Phys Med Rehabil 2012;91(8):689-96



- motor imagery, passive movement, movement observation, hand exercises in VR (Virtual Reality) can activate sensorimotor areas of the brain<sup>23,24,25,26,27</sup>.

Gloreha can be used in accordance with the most advanced rehabilitation treatments: Motor Imagery, Mirror Therapy, Constraint-induced Therapy, Bimanual Training.

---

<sup>22</sup> Sale P, Franceschini M. *Action observation and mirror neuron network: a tool for motor stroke rehabilitation*. Eur J Phys Rehabil Med. 2012 Jun;48(2):313-8. Epub 2012 Apr 20.

<sup>23</sup> Szameitat AJ, Shen S, Conforto A, Sterr A. *Cortical activation during executed, imagined, observed, and passive wrist movements in healthy volunteers and stroke patients*. Neuroimage. 2012 Aug 1;62(1):266-80. doi: 10.1016/j.neuroimage.2012.05.009. Epub 2012 May 11

<sup>24</sup> August K, Lewis JA, Chandar G, Merians A, Biswal B, Adamovich S. *fMRI analysis of neural mechanisms underlying rehabilitation in virtual reality: activating secondary motor areas*. Conf Proc IEEE Eng Med Biol Soc. 2006;1:3692-5.

<sup>25</sup> Tunik E, Saleh S, Adamovich SV. *Visuomotor discordance during visually-guided hand movement in virtual reality modulates sensorimotor cortical activity in healthy and hemiparetic subjects*. IEEE Trans Neural Syst Rehabil Eng. 2013 Mar;21(2):198-207. doi: 10.1109/TNSRE.2013.2238250. Epub 2013 Jan 9.

<sup>26</sup> Saleh S, Bagce H, Qiu Q, Fluet G, Merians A, Adamovich S, Tunik E. *Mechanisms of neural reorganization in chronic stroke subjects after virtual reality training*. Conf Proc IEEE Eng Med Biol Soc. 2011;2011:8118-21. doi: 10.1109/IEMBS.2011.6092002.

<sup>27</sup> Tunik E, Adamovich SV. *Remapping in the ipsilesional motor cortex after VR-based training: a pilot fMRI study*. Conf Proc IEEE Eng Med Biol Soc. 2009;2009:1139-42. doi: 10.1109/IEMBS.2009.5335392.